

# Gene Expression with the 192.24 IFC Using Delta Gene Assays

For more information, see the Real-Time PCR Analysis User Guide (PN 68000088) and the Juno System User Guide (PN 100-7070).

## Review Juno/IFC Controller RX Workflow

Load	Thermal-cycle (PCR) and image
Juno™ or RX	Biomark™ HD or Biomark

## Prepare the 192.24 IFC

### ! IMPORTANT


- Use the 192.24 Dynamic Array™ integrated fluidic circuit (IFC) within 24 hours of opening the package.
- Due to different accumulator volumes, use only syringes with 150 µL of control line fluid.
- Control line fluid on the IFC or in the inlets makes the IFC unusable.

- 1 Inject control line fluid into accumulator 2 (Acc2) on the IFC.
- 2 Remove and discard the blue protective film from the bottom of the IFC.

## Prepare 10X Assays

- 1 Prepare a 10X assay mix. Scale up appropriately for multiple runs.

Assay mix component	Vol. per Inlet (µL)	Vol. per inlet with overage (µL)	Vol. for 30 µL stock
100 µM each Delta Gene™ primers (forward and reverse combined), non-wet-lab tested (ASY-GE) or wet-lab tested (ASY-GE WET)	0.15	0.2	1.5
1X DNA suspension buffer	1.35	1.8	13.5
2X Assay Loading Reagent (PN 100-7611) ●	1.5	2.0	15.0
<b>Total</b>	<b>3.0</b>	<b>4.0</b>	<b>30.0</b>



## Prepare Sample Pre-Mix and Samples

- ! **IMPORTANT** Failure to do the following may result in a decrease in data quality.
- Pipet with care. The Delta Gene Sample Reagent is extremely viscous. **Do not vortex the Delta Gene Sample Reagent by itself at its stock concentration.**
  - Vortex thoroughly and centrifuge all assay and sample solutions **except** the 192.24 Delta Gene Sample Reagent before pipetting into IFC inlets. You can thaw the 192.24 Delta Gene Sample Reagent up to two times only.


- 1 Combine the SsoFast™ EvaGreen® Supermix with the 192.24 Delta Gene Sample Reagent. Vortex and spin down the sample pre-mix.
- 2 Aliquot 2.2 µL to each well of a 96-well plate. Next, add 1.8 µL of sample to individual wells. To remove any bubbles, vortex samples and centrifuge prior to adding the samples to the IFC.

Component	Vol. per inlet (µL)	Vol. per inlet with overage (µL)	Sample pre-mix with overage (µL)*
<b>SAMPLE PRE-MIX</b>			
2X SsoFast EvaGreen Supermix with low ROX™ (Bio-Rad PN 172-5211)	1.5	2.0	420.0
192.24 Delta Gene Sample Reagent (PN 100-6653) ●	0.15	0.2	42.0
Preamplified and Exo I-treated cDNA†	1.35	1.8	—
<b>Total</b>	<b>3.0</b>	<b>4.0</b>	<b>—</b>

\*For more information about PreAmp and Exonuclease I treatment, see Gene Expression PreAmp with Fluidigm PreAmp Master Mix and TaqMan Assays Quick Reference (PN 100-5875).

†For more information about PreAmp treatment, see Preamplification of cDNA for Gene Expression with Delta Gene Assays Quick Reference (PN 100-5875).

## 192.24 IFC Pipetting Map



## Load the IFC

### ! IMPORTANT

- While pipetting, do not go past the first stop on the pipette. Doing so may introduce air bubbles into inlets.
- For unused assay inlets, use 2.0  $\mu$ L assay loading reagent and 2.0  $\mu$ L of water per inlet.
- For unused sample inlets, use 2.2  $\mu$ L of sample pre-mix and 1.8  $\mu$ L of water per inlet.

- 1 Pipet 3  $\mu$ L of each assay and 3  $\mu$ L of each sample into the respective inlets on the IFC (see the 192.24 IFC pipetting map).
- 2 Pipet 150  $\mu$ L of Actuation Fluid (PN 100-6250) into the P1 well on the IFC.
- 3 Pipet 150  $\mu$ L of Pressure Fluid (PN 100-6249) into the P2 and P3 wells on the IFC.
- 4 Pipette 20  $\mu$ L of Pressure Fluid into the P4 and P5 wells on the IFC.

## For technical support visit [fluidigm.com/support](http://fluidigm.com/support)

EMAIL **United States** [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com) | **Europe** [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com)  
**Asia** [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com) | **Latin America** [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
**All other countries** [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com)

PHONE **United States (toll-free)** +1 866 358 4354 | **Europe** +33 160 92 42 40  
**Japan** +81 3 3662 2150 | **China (excluding Hong Kong)** +86 21 3255 8368  
**All other countries** +1 650 266 6100

## For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. Use standard laboratory safety protocols. **Safety data sheet information** [fluidigm.com/sds](http://fluidigm.com/sds) **Patent and license information** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices) **EU's WEEE directive information** [fluidigm.com/compliance](http://fluidigm.com/compliance) | Fluidigm, the Fluidigm logo, Biomark, Delta Gene, Dynamic Array, and Juno are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners.  
 © 2015 Fluidigm Corporation. All rights reserved. 06/2015

- 5 Place the IFC into the instrument and run the load script:
  - Juno: **Load Mix 192.24 GE**
  - RX: **Load Mix (169x)**

! **IMPORTANT** Start IFC run within 1 hour of loading samples.

## Collect Real-Time PCR Data

- 1 Remove any dust particles or debris from the IFC surface with clear tape.
- 1 Double-click the **Data Collection** icon on the desktop.
- 2 Select **Start a New Run**.
- 3 Verify that the camera and/or lamp are ready.
- 4 Confirm that the blue tape has been removed from the bottom of the IFC.
- 5 Select **Load**.
- 6 Verify IFC barcode and IFC type.
- 7 Choose project settings (if applicable). Click **Next**.
- 8 Provide a name and select a file storage location for a new IFC run, or browse to select a predefined run file. Click **Next**.
- 9 Choose the application, reference, and probes:
  - a Application type: **Gene Expression**
  - b Passive reference: **ROX**
  - c Assay: **Single probe**
  - d Probe type: **EvaGreen**
  - e Click **Next**.
- 10 Select thermal cycling protocol:
  - Biomark HD only (fast): **GE 192x24 Fast PCR+Melt v2.pcl**
  - Biomark HD or Biomark (standard): **GE 192x24 PCR+Melt v2.pcl**
- 11 Confirm **Auto Exposure** is selected.
- 12 Click **Next**.
- 13 Verify IFC run information.
- 14 Select **Start Run**.