

Real-Time PCR for Viral RNA Detection

This quick reference describes how to perform real-time PCR detection using the 192.24 Dynamic Array™ IFC (integrated fluidic circuit) with CDC assays (2019-nCoV CDC EUA Kit) from Integrated DNA Technologies (IDT™) on Biomark™ HD or Biomark. To prepare custom assays with this workflow, see the appendix on page 4.

IMPORTANT Before using this quick reference, read and understand the detailed instructions and safety guidelines in the Real-Time PCR for Viral RNA Detection Protocol (FLDM-00103).

Workflow

Workflow Step	Run Time*
1 Prepare and perform the RT reactions.	40 min
2 Preamplify the cDNA.	1 hr 40 min
3 Prepare the 192.24 IFC.	—
4 Prepare the CDC assays and sample mixes.	—
5 Load the IFC on Juno™ or IFC Controller RX.	30 min
6 Thermal-cycle and collect data on Biomark HD or Biomark.	Fast: 30 min Standard: 1 hr 30 min

* Does not include hands-on time

Prepare and Perform the RT Reactions

- 1 Thaw all reagents on ice. Briefly vortex and centrifuge the reagents before using.
- 2 On ice, prepare the reverse transcription (RT) reactions:

Table 1. RT reactions

Component	Vol per Reaction (μL)	Vol for 192 Reactions (μL)*
RNA template	5.0	—
Reverse Transcription Master Mix (100-6297)	1.25	264
Total	6.25	—

* Includes overage

- a Add 5.0 μL of RNA template into each well of 2 new 96-well PCR plates (on ice).
 - b Aliquot 33 μL of Reverse Transcription Master Mix into each well of an 8-well strip.
 - c Use an 8-channel pipette to transfer 1.25 μL of Reverse Transcription Master Mix into each well containing RNA in the 96-well plates from Step a.
- 3 Properly seal and gently vortex to mix the RT reactions.

- 4 Centrifuge the reactions at 3,000 × g for 60 sec, then place the plates in a standard thermal cycler for 96-well plates and incubate using the following protocol:

Temperature	Time	Condition
25 °C	5 min	Hold
42 °C	30 min	Hold
85 °C	5 min	Hold
4 °C	∞	Hold

STOPPING POINT After the RT reaction is complete, the reactions can be stored at -20 °C or used immediately for preamplification reactions with Preamp Master Mix.

Preamplify the cDNA

NOTE If you are using custom assays, see the appendix on page 4 for the assay preparation and pooling procedure.

Pool the CDC Assays

- 1 Briefly vortex and centrifuge the reagents before use.
- 2 Dilute the pooled CDC assays as shown in Table 2.

Table 2. Pooled CDC assays

Component	Vol for 3 Assays/ 192 Reactions (μL)
CDC assays (6.7 μM primer, forward and reverse; 1.7 μM probe) (N1, N2, and RNase P)	10.5 μL × each assay = 31.5
Dilution Reagent (100-8730)	668.5
Total	700

Final concentration of pooled CDC assays:
100.5 nM primer; 25.5 nM probe

NOTE Volume can be adjusted proportionally based on the number of samples to be amplified.

Prepare the Preamplification Sample Mix

- 1 Briefly vortex and centrifuge the reagents before use.
- 2 In a DNA-free hood, combine the components (Table 3) in a new 1.5 mL tube to make the preamplification pre-mix. Scale up appropriately for multiple runs.

Table 3. Preamplification pre-mix

Component	Vol per Reaction (μL)	Vol for 192 Reactions/ 1 IFC (μL)*
Preamp Master Mix (100-5744) ●	2.5	528
Pooled CDC assays (see Table 2)* (100.5 nM primer, forward and reverse; 25.5 nM probe)	3.125	660
PCR Water (100-5941)	0.625	132
Total	6.25	1,320

Final concentration of CDC assays in preamplification pre-mix: 25.125 nM primer; 6.375 nM probe

* If using custom assays, see Table 9 on page 4.

† Includes overage

- 3 Cap tube, vortex and centrifuge the preamplification pre-mix.
- 4 Add preamplification pre-mix to the RT reactions to prepare the preamplification sample mix (Table 4).

Table 4. Preamplification sample mix

Component	Vol per Inlet (μL)
Preamplification sample pre-mix (see Table 3)	6.25
RT reactions (see Prepare and Perform the RT Reactions on page 1)	6.25
Total	12.5

- a Aliquot 165 μL of preamplification sample pre-mix into each well of an 8-well strip.
- b Use an 8-channel pipette to transfer 6.25 μL of pre-mix into each well containing the RT reactions from Prepare and Perform the RT Reactions on page 1.
- 5 Seal the plates, then vortex and centrifuge them at 3,000 × g for 60 sec.

Thermal-Cycle the Preamplification Sample Mix

- 1 Place the plates in a standard thermal cycler for 96-well plates and cycle using the following table as a guide:

Temperature	Time	Condition	
95 °C	2 min	Hold	
95 °C	15 sec	20 cycles	Denaturation
60 °C	4 min		Annealing/extension
4 °C	∞	Hold	

NOTE The appropriate number of cycles may need to be determined empirically.

- 2 After cycling, dilute the reaction 1:5 by adding 50 μL of Dilution Reagent (100-8730) to the final 12.5 μL reaction volume for a total volume of 62.5 μL.

STOPPING POINT Diluted reaction products (preamplified cDNA) can either be assayed immediately or stored at -20 °C for later use. Diluted reaction products should be stable for at least 1 week.

Prepare the 192.24 IFC

IMPORTANT

- Use the IFC within 24 hr of opening the package.
- Due to different accumulator volumes, only use 192.24 syringes with 150 μL of Control Line Fluid (100-4058).
- Control Line Fluid on the IFC or in the inlets makes IFC unusable.

- 1 Ensure that the RX Interface Plate is installed in the Juno instrument.
- 2 Inject Control Line Fluid into the top accumulator on the IFC (Figure 1). Use the entire contents of the syringe.
- 3 Remove and discard the protective film from bottom of the IFC.




Figure 1. Top accumulator on 192.24 IFC

Prepare the CDC Assay Mixes

NOTE If you are using custom assays, see the appendix on page 4 for the assay mixes preparation procedure.

- 1 Briefly vortex and centrifuge the reagents before use.
- 2 In a DNA-free hood, prepare aliquots of 10X CDC assay mixes using components in Table 5. Scale up appropriately for multiple runs.

Table 5. CDC assay mixes

Component	Vol per Inlet (μL)†	Vol for N Assays (μL)
CDC assays (2019-nCoV CDC EUA Kit)* (IDT, 10006606)	3.0	3.0 × N
4X Assay Loading Reagent (102-0114) ●	1.0	1.0 × N
Total	4.0	4.0 × N

Final concentration (at 10X): 5 μM primers; 1.25 μM probe

* For unused assay inlets, replace the CDC assays with 3.0 μL of PCR Water (100-5941).

† Includes overage

Prepare the Sample Mix

- 1 Briefly vortex and centrifuge the reagents before use.
- 2 In a DNA-free hood, combine the components (Table 6) in a sterile 1.5 mL tube to make the sample pre-mix. Scale up appropriately for multiple runs.

NOTE This is enough volume to fill the entire IFC.

Table 6. Sample pre-mix

Component	Vol per Inlet (μL)*	Sample Pre-Mix for 192.24 IFC (μL)†
Master Mix		
• For standard: TaqMan Gene Expression Master Mix (Thermo Fisher Scientific, 4369016)	2.0	440.0
• For fast: TaqMan Fast Advanced Master Mix (Thermo Fisher Scientific, 4444557)		
20X GE Sample Loading Reagent (85000735)	0.2	44.0
Total	2.2	484.0

* Includes overage

† 220 reactions for ease of pipetting

a Aliquot 60.5 μL of sample pre-mix into each well of a new 8-well strip

b Use an 8-channel pipette to transfer 2.2 μL of pre-mix into each well of 2 new 96-well plates.

- 3 Remove the plates from the DNA-free hood and prepare the final sample mix by adding preamplified cDNA to each well (Table 7).

Table 7. Final sample mix

Component	Vol per Inlet (μL)†
Sample pre-mix (see Table 6)	2.2
Preamplified cDNA* (see Preamplify the cDNA on page 1)	1.8
Total	4.0

* For unused sample inlets, replace the preamplified cDNA with 1.8 μL of PCR Water (100-5941).

† Includes overage

- 4 Seal the plates, then vortex and centrifuge them at 3,000 × g for 60 sec.

Load the IFC

IMPORTANT

- Vortex thoroughly and centrifuge all assay and sample solutions before pipetting into the IFC inlets. Failure to do so may result in a decrease in data quality.
- While pipetting, do not go past the first stop on the pipette. Doing so may introduce air bubbles into inlets.

- 1 Pipet 3 μL of each sample into the respective inlets on the IFC.
- 2 Pipet 3 μL of each assay into the respective inlets on the IFC.
- 3 Pipet 150 μL of Actuation Fluid (100-6250) into the P1 reservoir () on the IFC.
- 4 Pipet 150 μL of Pressure Fluid (100-6249) into each of the P2 and P3 reservoirs () on the IFC.
- 5 Pipet 20 μL of Pressure Fluid into each of the P4 and P5 inlets (○) on the IFC.
- 6 Blot the IFC surface with a dry, lint-free cloth.

- 7 Place the IFC into the controller:

- Juno: Tap **OPEN** to open the instrument tray and align the notched corner of the IFC to the white notch on the tray. Tap **LOAD**.
- RX: Press **EJECT** to open the instrument tray and align the notched corner of the IFC to the A1 mark. Press **Load Chip**.

- 8 Run the Load Mix script:

- Juno: Tap **Load Mix 192.24 GE**, then tap **Run**.
- RX: Select **Load Mix (169x)** and press **Run Script**.

IMPORTANT Start the IFC run within 1 hr of loading samples.

- 9 If necessary, turn on the Biomark system (computer and instrument), launch the Data Collection software, and turn on the lamp. The lamp takes 20 min to warm up.




Figure 2. Pipetting map for the 192.24 IFC

Collect Data

Transfer the loaded IFC from the controller to the Biomark HD or Biomark (standard only) instrument and collect data.

- 1 Use clear tape to remove any dust particles or debris from the IFC surface.
- 2 If necessary, double-click the **Data Collection** icon on the desktop of the Biomark HD or Biomark system computer to launch the software.
- 3 Click **Start a New Run**.
- 4 Confirm that the camera status indicator and the lamp status indicator (Biomark only) at the bottom of the window are green.
- 5 Place the loaded IFC on the instrument tray and align the notched A1 corner on the IFC with the A1 label on the tray. Click **Load**.
- 6 Verify the IFC barcode and IFC type and then click **Next**.
- 7 Complete the Chip Run section by selecting either a new or a pre-defined run.

- 8 Complete the Chip Run Name and Location section:
 - a Enter a run name or select the checkbox to use the IFC barcode as the run name.
 - b Select a file storage location for a new IFC run or browse to select a pre-defined run file and click **Next**.

- 9 Complete the Application, Reference and Probes section and then click **Next**.

For	Select
Application	Gene Expression
Passive reference	ROX™
Assay	Single probe
Probes	FAM-MGB

- 10 Browse to and select the thermal protocol:

For	Select
Standard	GE 192x24 Standard v1.pcl
Fast (Biomark HD only)	GE 192x24 Fast v1.pcl

- 11 Confirm that **Auto Exposure** is selected. Click **Next**.
- 12 Verify the IFC run information and click **Start Run**.
- 13 After the run is complete, analyze your data using the Real-Time PCR Analysis software.

Appendix: Custom Assays

Prepare the 20X Custom Assays

Based on the final reaction concentration of primers and probes, prepare the 20X custom assays. For example, Table 8 shows the preparation for a reaction with a final target of concentration for forward and reverse primers each at 500 nM and probe at 125 nM.

Table 8. 20X custom assays

Component	Vol per Inlet (μL)	Final 20X Concentration (μM)
Forward primer (100 μM)	10	10
Reverse primer (100 μM)	10	10
TaqMan probe (FAM-MGB) (100 μM)	2.5	2.5
Dilution Reagent (100-8730) <input type="radio"/>	77.5	—
Total	100	—

Pool the 20X Custom Assays for Preamplication

- 1 Briefly vortex and centrifuge the reagents before use.
- 2 In a microcentrifuge tube, combine equal volumes of each 20X custom assay for preamplication, up to a total of 24 assays.

- 3 Dilute the pooled custom assays as shown in Table 9. Table 9 provides an example using 24 assays.

NOTE The primer concentration in the preamplication reaction may need optimization.

Table 9. Pooled custom assays for preamplication

Component		Vol for 24 Assays/192 Reactions (μL)*
20X assays	7 μL × each assay	168
Dilution Reagent (100-8730) <input type="radio"/>	700 – (7 μL × each assay)	532
Total		700

* Includes overage

- 4 Go to [Prepare the Preamplication Sample Mix on page 1](#).

Prepare the 10X Custom Assay Mixes

The 10X custom assay mixes require 2X Assay Loading Reagent, included in the 192.24 GE Dynamic Array Reagent Kit (100-6267). If you are using the 192.24 GE Dynamic Array 4X Reagent Kit (102-0166), obtain 2X Assay Loading Reagent by either:

- Ordering from fluidigm.com (100-7611), or
- Preparing by combining components (Table 10) in a sterile tube.

Table 10. 2X Assay Loading Reagent

Component	Vol (μL)
4X Assay Loading Reagent (102-0114)	30
PCR Water (100-5941)	30
Total	60

- 1 Briefly vortex and centrifuge reagents before use.
- 2 In a DNA-free hood, prepare aliquots of 10X custom assay mixes using the components in Table 11. Scale up appropriately for multiple runs.

Table 11. 10X custom assay mixes

Component	Vol per Inlet (μL) [†]	Vol for 24 Assays (μL)
20X custom assays (see Table 8)*	2.0	48.0
2X Assay Loading Reagent	2.0	48.0
Total	4.0	96.0

* For unused assay inlets, replace the 20X assays with 2.0 μL of PCR Water (100-5941).

[†] Includes overage

NOTE Volume can be adjusted proportionally based on the number of samples to be amplified.

- 3 Go to [Prepare the Sample Mix on page 2](#).

For technical support visit techsupport.fluidigm.com.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com
 Europe/Middle East/Africa/Russia +44 1223 598100 | techsupport@fluidigm.com
 Japan +81 3 3662 2150 | techsupport@fluidigm.com

Latin America +1 650 266 6100 | techsupport@fluidigm.com

China (excluding Hong Kong) +86 21 3255 8368 | techsupport@fluidigm.com

All other Asian countries/India/Australia +1 650 266 6100 | techsupport@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legal/notices. **Limited Use License to Perform Preamplication with Fluidigm IFCs:** A license to use Thermo Fisher Scientific's patented preamplication method workflows involving a Fluidigm integrated fluidic circuit (IFC) can be obtained (i) with purchase of a Fluidigm IFC from Fluidigm Corporation or (ii) by a separate license from Thermo Fisher Scientific. For licensing information, contact outlicensing@lifetech.com. **Trademarks:** Fluidigm, the Fluidigm logo, Biomark, Dynamic Array, EP1, and Juno are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2020 Fluidigm Corporation. All rights reserved. 04/2020