Safety Information for Mass Cytometry Systems

For use with CyTOF® 2, Helios™, and the Hyperion™ Imaging System

Contents

Safety Guide	2
Introduction	2
Safety Guidelines	3
Sikkerhedsvejledning	10
Indledning	10
Sikkerhedsretningslinjer	11
Sicherheitsleitfaden	18
Einleitung	18
Sicherheitsleitlinien	19
Guía de seguridad	27
Introducción	27
Pautas de seguridad	28
Guide de sécurité	36
Introduction	36
Directives de sécurité	37
Guida alla sicurezza	45
Introduzione	45
Linee guida sulla sicurezza	46
Veiligheidshandleiding	54
Inleiding	54
Veiligheidsvoorschriften	55

Safety Guide

For Research Use Only. Not for use in diagnostic procedures.

Introduction

This document provides important safety information pertaining to the operation of Fluidigm® mass cytometry systems. For detailed information about using a system, see the appropriate user guide.

IMPORTANT Before using the system, read and understand the safety guidelines in this document. Failure to follow these guidelines may result in undesirable effects, injury to personnel, and/or damage to the system or to property.

Operator Safety for the Hyperion Imaging System

The Hyperion™ Imaging System is classified as a Class 1 laser device.

The laser radiation of the Class 1 laser system is eye-safe under all operating conditions. This product includes lasers of a higher class whose beams are confined within a suitable enclosure so that access to laser radiation is physically prevented.

Safety Alert Conventions

Fluidigm documentation uses specific conventions for presenting information that may require your attention. Refer to the following safety alert conventions.

Safety Alerts for Chemicals

For hazards associated with chemicals, this document follows the United Nations Globally Harmonized System of Classification and Labelling of Chemicals (GHS) and uses indicators that include a pictogram and a signal word that indicates the severity level:

Indicator	Description
	Pictogram (see example) consisting of a symbol on a white background within a red diamond-shaped frame. Refer to the individual safety data sheet (SDS) for the applicable pictograms and hazards pertaining to the chemicals being used.
DANGER	Signal word that indicates more severe hazards.
WARNING	Signal word that indicates less severe hazards.

Safety Alerts for Instruments

For hazards associated with instruments, this document uses indicators that include a pictogram and signal words that indicate the severity level:

Indicator	Description
	Pictogram (see example) consisting of a symbol on a white background within a black triangle-shaped frame. Refer to the system user guide for the applicable pictograms and hazards pertaining to system usage.
DANGER	Signal word that indicates an imminent hazard that will result in severe injury or death if not avoided.
WARNING	Signal word that indicates a potentially hazardous situation that could result in serious injury or death if not avoided.
CAUTION	Signal word that indicates a potentially hazardous situation that could result in minor or moderate personal injury if not avoided.
IMPORTANT	Signal word that indicates information necessary for proper use of products or successful outcome of experiments.

Safety Data Sheets

Read and understand the SDSs before handling chemicals. To obtain SDSs for chemicals ordered from Fluidigm, either alone or as part of this system, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

Some chemicals referred to in this user guide may not have been provided with your system. Obtain the SDSs for chemicals provided by other manufacturers from those manufacturers.

Safety Guidelines

General Safety

In addition to your site-specific safety requirements, Fluidigm recommends the following general safety guidelines in all laboratory and manufacturing areas:

- Use the appropriate personal protective equipment (PPE): safety glasses, fully enclosed shoes, lab coats, and gloves, according to your laboratory safety practices.
- Know the locations of all safety equipment (fire extinguishers, spill kits, eyewashes/showers, first-aid kits, safety data sheets, etc.), emergency exit locations, and emergency/injury reporting procedures.
- Do not eat, drink, or smoke in lab areas.
- Maintain clean work areas.
- Wash hands before leaving the lab.

Laser Safety for the Hyperion Tissue Imager

The Hyperion Tissue Imager is a solid-state pulsed laser classified as a Class 1 laser system. It complies with 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007.

Instrument Safety

The system should be serviced by authorized personnel only.

WARNING Do not modify this instrument or system. Unauthorized modifications may create a safety hazard.

WARNING BIOHAZARD. If you are putting biohazardous material on the instrument or system, use appropriate personal protective equipment and adhere to Biosafety in Microbiological and Biomedical Laboratories (BMBL), a publication from the Centers for Disease Control and Prevention, and to your lab's safety protocol to limit biohazard risks. If biohazardous materials are used, properly label the equipment as a biohazard. For more information, see the BMBL guidelines online at cdc.gov/biosafety/publications/index.htm.

WARNING Laser instruments generate potentially hazardous UV radiation. Do not remove the top cover of the Hyperion Tissue Imager. Only a Fluidigm field service engineer should remove the top cover and perform maintenance.

WARNING TIPPING HAZARD. The Hyperion Tissue Imager has a high center of gravity, and therefore there is risk of tipping when moving the instrument or system.

WARNING HOT SURFACE HAZARD. A safety interlock on the CyTOF 2 and Helios systems automatically shuts off the plasma if the chamber and interface are not fully coupled. Do not defeat the interlock. Do not remove the shield that protects the sample introduction system. The heat shield is designed to protect users from burns from the heater.

WARNING HOT SURFACE HAZARD. The torch components, the vacuum interface, and the sample introduction system components remain hot for some time after the plasma has been shut off. Allow sufficient time for these items to cool to room temperature before you handle them.

WARNING PHYSICAL INJURY HAZARD. When installing or moving the instrument or system, contact a Fluidigm field service engineer for assistance. See the user guide for the weight of the boxed or crated instrument or system.

CAUTION FINGER CUT HAZARD. The injector may break if excessive pressure is used to remove the injector from the injector holder.

CAUTION FINGER CUT HAZARD. Broken glass may cause injury or cutting of fingers or hands. Use caution when loading and unloading the sample slides and when handling instrument glassware, including the injector and torch holder/body.

Symbols on the Instrument

The following table describes the hazard symbols that may be used in this document or on labels on the system.

Symbol	Description
	Hazard. Consult the user guide for further information.
	Hot surface hazard. Do not touch; potential for personal injury.

Symbol	Description
	Biohazard.
	Electricity hazard. Indicates high electricity levels and a threat of electric shock from machines and/or equipment in the vicinity. You may suffer severe injuries or death.
	Pinch hazard. Indicates where pinch hazards exist. Exercise caution when operating around these areas.
	
	Lifting hazard.
	Indicates rotating blades can crush or cut fingers or hands. Keep hands clear.
	Laser hazard. Indicates the presence of a laser.
	Finger cut hazard. Broken glass may cause injury or cutting of fingers or hands. Caution when loading and unloading the sample slides.
	Non-ionizing radiation hazard. Exposure to high-frequency radio waves and radio frequency radiation can result in injuries.
	Tipping hazard. Movement or impact with the instrument or system may cause tipping.
	Trip hazard. Watch your step to avoid falling over objects.
	Indicates specific chemical harm.
	Indicates hazardous, toxic, or very toxic materials that are very hazardous to health or potentially fatal when inhaled, swallowed, or in contact with the skin.
	Indicates caustic and acid materials that can destroy the skin and eat through metals.
	Indicates the presence of material contained under pressure, including compressed gas, dissolved gas, or gas liquefied by compression or refrigeration.
	A compressed gas cylinder can become a projectile when ruptured, with the potential to cause significant damage.
	Indicates a health hazard.

Symbol	Description
	Power and standby symbol.
	Power switch is in the Off position.
	Power switch is in the On position.
	Protective conductor terminal (main ground). It must be connected to earth ground before any other electrical connections are made to the instrument or system.
	To minimize negative environmental impact from disposal of electronic waste, do not dispose of electronic waste in unsorted municipal waste. Follow local municipal waste ordinances for proper disposal provision. Contact customer service for information about responsible disposal options.

Electrical Safety

WARNING ELECTRICAL HAZARD. DO NOT REMOVE THE COVERS. Electrical shock can result if the system is operated without its protective covers. No internal components under the covers are serviceable by the user.

WARNING ELECTRICAL HAZARD. Plug the system into a properly grounded receptacle with adequate current capacity.

WARNING Lethal voltages are present at certain areas within the system. Only a Fluidigm field service engineer or those similarly authorized and trained by Fluidigm personnel should install or repair the system.

WARNING Do not remove the side panel on the electrical box of the Hyperion Tissue Imager. Only a Fluidigm field service engineer should remove the side panel and maintain the electrical box.

WARNING The system power supplies are capable of generating potentially lethal voltages and currents. Store the removable system handle separately from the system. Maintenance should be performed only by a Fluidigm field service engineer or by maintenance personnel, employed by the customer, who have been trained by Fluidigm and are appropriately certified.

WARNING The interface and backing pumps in the system are in close proximity to areas where high voltages are present. User access to the pumps is not advised. Only Helios™ operators trained by Fluidigm may access the pump compartment and change the interface pump oil. Disengage the radio frequency generator circuit breaker on the right side of the system before accessing this area.

WARNING Do not touch electrical wires, contacts, transformers, or transformer components during the oil inspection procedure (see the user guide). A metal shield located in the system compartment above the interface pump contains the transformers and electrical wires. There is no need to access this section when servicing the pumps.

- When the system is connected to line power, opening system covers is likely to expose live electrical parts.
- High voltage can still be present even when the radio frequency generator power switch is off.
- Water lines should be located away from electrical connections. Condensation and potential leaks may create an unsafe environment in the proximity of electrical connections.

WARNING Before performing maintenance on the cones or torch, switch off the radio frequency generator power using the circuit breaker at the right rear of the system. Wait at least 5 minutes for residual electrical charge to dissipate. Additional time is required to allow the inductively coupled plasma torch, cones, and load coil to reach room temperature.

- Capacitors inside the system may still be charged even if the system has been disconnected from all voltage sources.
- The system must be correctly connected to a suitable electrical supply (see the site requirements guide for further details).
- The power supply must have a correctly installed protective conductor (earth ground) and must be installed or checked by a qualified electrician before connecting the system.

WARNING Any interruption of the protective conductor (earth ground) inside or outside the system or disconnection of the protective conductor terminal is likely to make the system dangerous.

- Do not operate the system with any covers or internal parts removed.
- Do not attempt to perform internal adjustments or replacements except as directed in this user guide.

Chemical Safety

The responsible individuals must take the necessary precautions to ensure that the surrounding workplace is safe and that system operators are not exposed to hazardous levels of toxic substances. When working with any chemicals, refer to the applicable safety data sheets (SDSs) provided by the manufacturer or supplier. When handling any chemical, the following safe-handling guidelines should be strictly observed:

- Do not inhale fumes from chemicals. Use adequate ventilation and return caps to bottles immediately after use.
- Use, store, and dispose of chemicals according to manufacturer recommendations and to regulations applicable to the locality, state, province, and/or country.
- When preparing chemical solutions, always work in a fume hood that is suitable for those chemicals.
- Conduct sample preparation away from the system to minimize corrosion and contamination.
- Store solvents in an approved cabinet (with the appropriate ventilation) away from the system.

Laboratory Ventilation Safety

Toxic combustion products, metal vapor, and ozone can be generated by the system, depending on the type of analysis. An efficient ventilation system must be provided for your system. When the plasma is on, hot gases are vented through two exhaust vents at the back of the system. Detailed information on exhaust vents is included in the site requirements guide.

WARNING Use of the instrument or system without adequate ventilation to outside air may constitute a health hazard. Take extreme care to vent exhaust gases properly.

WARNING The instrument or system is designed for analysis of fixed/permeabilized, non-live cells only. Under normal operation, cells are completely combusted in the inductively coupled plasma. High levels of UV radiation inside the torch box are significantly above lethal levels for most single airborne cells. However, in the event of plasma shutdown, the non-ionized portion of a sample can enter the torch box exhaust gases.

Drain Vessel Safety for CyTOF 2

A drain vessel is supplied with the CyTOF® 2 instrument. The vessel is made of high-density polyethylene (HDPE) and is used to gather the effluent from the Flow Injection Valve of the sample introduction system. For safe operation of your system, properly install and maintain the drain vessel and drain tubing. Waste disposal procedures must be in accordance with all national, state/provincial, and local health and safety regulations and laws. Drain vessels may contain flammable, acidic, caustic, or organic solutions; cell debris; and small amounts of the elements analyzed.

WARNING It is necessary to follow appropriate waste segregation guidelines in order to prevent effluents from reacting in the drain vessel.

- Never place the vessel in an enclosed cabinet. Doing so may result in a buildup of hazardous gases.
- Do not use a glass drain vessel. A glass drain vessel may break and spill toxic or corrosive liquids.
- Place the drain vessel in an area that is visible to the operator, who can observe the level of collected effluent and empty the vessel when necessary.
- Check the condition of the drain tubing regularly to monitor deterioration. Organic solvents deteriorate the tubing more quickly than aqueous solutions. When the tubing becomes brittle or cracked, replace it.
- Empty the drain bottle regularly. Dispose of waste in accordance with all national, state/provincial, and local health and safety regulations and laws.

Pressurized Gas Safety

Safe Handling of Gas Cylinders

Argon gas used with the system is normally stored in liquid argon tanks or pressurized containers. Carefully use, store, and handle compressed gases in cylinders. Gas cylinders can be hazardous if they are mishandled. Argon is neither explosive nor combustible.

Helium gas is supplied in the non-liquefied or liquid form in a compressed gas cylinder for use with the Hyperion™ Imaging System.

Contact the gas supplier for a safety data sheet containing detailed information on the potential hazards associated with the gas.

IMPORTANT If liquid argon or liquid helium is used, the gas cylinder must be fitted with an overpressure regulator, which will vent the cylinder as necessary to prevent it from becoming a safety hazard.

WARNING Do not use electronic pressure regulator and auto switching valves because doing so may affect the plasma stability and may result in frequent loss of plasma.

WARNING It is recommended to install an oxygen sensor in the room where the operator and gas storage are located.

Sample Handling and Preparation Safety

WARNING For better control of contamination, dedicate laboratory reagents and consumables to use with CyTOF® instruments and Maxpar® reagents only.

Radio Frequency Radiation Safety

WARNING RADIO FREQUENCY RADIATION. The system generates high levels of RF energy, which is potentially hazardous if allowed to escape. The system is designed to contain the RF energy within the shielded enclosures of the torch compartment and the RF power supply. Safety interlocks prevent the system from operating without all covers, doors, and shields in place.

WARNING Inductively coupled plasma-based systems generate high levels of radio frequency (RF) energy within the RF power supply and the torch box. RF energy is potentially hazardous if allowed to escape. Do not bypass or disconnect safety devices and safety interlocks.

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For technical support visit techsupport.fluidigm.com.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | support.northamerica@fluidigm.com Latin America +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 598100 | support.europe@fluidigm.com China (excluding Hong Kong) +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com All other Asian countries/India/Australia +1 650 266 6100 | techsupportasia@fluidigm.com

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Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legal/notices. **Limited Use Label License:** The purchase of this Fluidigm Instrument and/or Consumable product conveys to the purchaser the limited, nontransferable right to use with only Fluidigm Consumables and/or Instruments respectively except as approved in writing by Fluidigm. **Trademarks:** Fluidigm, the Fluidigm logo, CyTOF, Helios, Hyperion, and Maxpar are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2019 Fluidigm Corporation. All rights reserved. 11/2019 GEN-SG-1