

Purchases

Beckman Coulter

RNAClean XP beads

A63987

AMPure XP Beads

A63880

Life Technologies

Arraysript (AM2048 or from MessageAmplII kit – AM1751)

NEB

Phusion[®] High-Fidelity PCR Master Mix with HF Buffer – M0531

Primer manufacturer

randomhexRT primer – GCCTTGGCACCCGAGAATTCCANNNNNN – used at 100 μ M

small RNA PCR primers - PCR primers used at 10 μ M.

RP1: 5'-AATGATACGGCGACCACCGAGATCTACACGTTTCAGAGTTCTACAGTCCG-s-A-3'

RPI1: 5'-CAAGCAGAAGACGGGCATACGAGATCGTGATGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI2: 5'-CAAGCAGAAGACGGGCATACGAGATACATCGGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI3: 5'-CAAGCAGAAGACGGGCATACGAGATGCCTAAGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI4: 5'-CAAGCAGAAGACGGGCATACGAGATTGGTCAGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI5: 5'-CAAGCAGAAGACGGGCATACGAGATCACTGTGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI6: 5'-CAAGCAGAAGACGGGCATACGAGATATTGGCGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI7: 5'-CAAGCAGAAGACGGGCATACGAGATGATCTGGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI8: 5'-CAAGCAGAAGACGGGCATACGAGATCAAGTGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI9: 5'-CAAGCAGAAGACGGGCATACGAGATCTGATCGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI10: 5'-CAAGCAGAAGACGGGCATACGAGATAAGCTAGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI11: 5'-CAAGCAGAAGACGGGCATACGAGATGTAGCCGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

RPI12: 5'-CAAGCAGAAGACGGGCATACGAGATTACAAGGTGACTGGAGTTCCTTGGCACCCGAGAATTCC-s-A-3'

aRNA Clean up and fragmentation

Clean up

- 1) Pool 5 μ l of each of the harvested aRNA libraries (~480 μ l total)
- 2) Split into four 0.5-ml tubes (~120 μ l each)
- 3) Add 1.5x volumes (180 μ l) of RNAClean XP beads to each. Mix well until the liquid appears homogeneous. Incubate for 5 min at room temperature.
- 4) Put the tubes on a magnetic stand for bead separation for 5 min.
- 5) Aspirate and discard the supernatant.
- 6) Wash 3 times with 70% ethanol.
- 7) Remove the ethanol and let dry for 10 min or until completely dry.
- 8) Elute the aRNA with 20 μ l TE or RNase-free DDW. Combine eluates and mix until homogeneous.
- 9) Place the tube on a magnetic stand for 5 min
- 10) Transfer the eluted aRNA to a new Eppendorf

Fragmentation

To 20µl aRNA add 5µl of fragmentation buffer (200mM Tris-acetate, pH 8.1, 500 mM KOAc, 150 mM MgOAc)

Incubate for 3 min at 94°C

Move to ice and add 2.5µl stop solution (0.5 M EDTA pH8)

Repeat the clean up as above, but add 1× beads (27.5µl). Elute in 20µl.

Check aRNA amount and quality on 2100 Bioanalyzer

Library preparation

Protocol designed for 5-10ng amplified RNA, although as little as 1-2ng can be used, but then additional PCR cycles are required. Sample volume should be adjusted to 5µl, either by adding water or drying down in a speedvac, depending on RNA concentration.

RT reaction (Ambion kit)

To 5 µl aRNA add 1 µl randomhexRT primer.

Incubate 2min at 70°C, quick chill on ice.

Add 4µl of the following mix to each reaction:

First Strand buffer	1µl
dNTP (2.5mM)	2µl
RNase Inhibitor	0.5µl
ArrayScript	0.5µl

Incubate 10 min at 25°C.

Incubate 2hr at 42°C (in hybridization oven, or thermal cycler with lid at 50°C)

PCR amplification:

To each reverse transcription reaction add 38µl of the following mix:

- Ultra Pure Water 11µL
- Phusion PCR mix 25µL
- RNA PCR Primer (RP1, from Illumina small RNA kit) 2µL

To each reaction add 2µl of a uniquely indexed RNA PCR Primer (RPIX, from Illumina small RNA kit)

Amplify the tube in the thermal cycler using the following PCR cycling conditions:

30 seconds at 98°C

12 cycles of:

10 seconds at 98°C

30 seconds at 60°C

30 seconds at 72°C

10 minutes at 72°C

Hold at 4°C

Can go up to 15 cycles if necessary, or down to 11 if starting with the full 10ng.

Bead Cleanup of PCR products – Repeat 1:

- Prewarm beads to room temperature.
- Vortex AMPure XP Beads until well dispersed, then add 50 μ l to the 50 μ l PCR reaction. Mix entire volume up ten times to mix thoroughly.
- Incubate at room temperature for 15 min.
- Place on magnetic stand for at least 5 min, until liquid appears clear.
- Remove and discard 95 μ l of the supernatant.
- Add 200 μ l freshly prepared 80% EtOH.
- Incubate at least 30 seconds, then remove and discard supernatant without disturbing beads.
- Add 200 μ l freshly prepared 80% EtOH
- Incubate at least 30 seconds, then remove and discard supernatant without disturbing beads.
- Air dry beads for 15 min, or until completely dry.
- Resuspend with 30 μ l Resuspension Buffer (from Illumina kit). Pipette entire volume up and down ten times to mix thoroughly.
- Incubate at room temperature for 2 min.
- Place on magnetic stand for 5 min, until liquid appears clear.
- Transfer 30 μ l of supernatant to new tube.

Bead Cleanup of PCR products – Repeat 2:

Repeat as above, but adding 24 μ l beads and elute in 10 μ l Resuspension Buffer at the end, transferring 10 μ l to a new tube.