

Juno System Software

To download the latest version of the software, go to fluidigm.com/software.

v3.14.1

v3.14.1 Enhancement

48.Atlas™ IFC support. New scripts to run the 48.Atlas IFC (integrated fluidic circuit).

v3.14.1 Bug Fixes

One Step LP-8.8.6 script. The required thermal protocols are now included when exporting the script to a USB drive.

v3.13.1

v3.13.1 Enhancements

LP 8.8.6 IFC support. New scripts to run the LP 8.8.6 IFC.

v3.13.1 Bug Fixes

None.

v3.12.1

v3.12.1 Enhancements

- Added compatibility with the Microsoft® Windows® 10 operating system.
- Moved the instrument log folder to C:\Fluidigm\InstrumentLog.
- The About This System page displays the software build number along with the software version.

v3.12.1 Bug Fixes

None.

v3.11.1

v3.11.1 Enhancements

- **24.192 GE IFC support.** New scripts to run the 24.192 Dynamic Array™ IFC for Gene Expression.
- **192.24 GE and 192.24 GT IFC support.** Added Prime to the description of the Load Mix scripts for the 192.24 Gene Expression and 192.24 Genotyping IFCs.

v3.11.1 Bug Fixes

None.

v3.9.1

v3.9.1 Enhancements

None.

v3.9.1 Bug Fixes

Switch interface plate. Correct display of TX Interface Plate type when switching the interface plate.

v3.7.1

v3.7.1 Enhancements

None.

v3.7.1 Bug Fixes

IFC orientation check. Orientation for IFCs is now determined by the location of the barcode label. If the barcode is not detected, the user is prompted to check the orientation of the IFC.

v3.5.1

v3.5.1 Enhancement

LP 48.48 IFC support. New scripts to run the LP 48.48 IFC.

v3.5.1 Bug Fixes

- **Leap year.** Scripts recognize 29 February as a valid date.
- **Daylight saving time.** Instrument clock adjusts for daylight saving time.

v3.3.4

v3.3.4 Enhancement

LP 192.24 IFC support. New scripts to run the LP 192.24 IFC.

v3.3.4 Bug Fixes

None.

v3.1.4

v3.1.4 Enhancements

- **System power down.** This new operation on Juno™ permits safe shutdown of the operating system before powering down the instrument. While the system is running, press the front power and standby indicator for 2 seconds to initiate the shutdown.
- **Barcode recognition.** Improved to recognize all supported IFCs on the instrument.
- **System cleaning script.** Updated with a running time of 13 minutes.

v3.1.4 Bug Fixes

None.

For technical support visit techsupport.fluidigm.com.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com Latin America +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 598100 | techsupporteurope@fluidigm.com China (excluding Hong Kong) +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com All other Asian countries/India/Australia +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legalnotices. **Trademarks:** Fluidigm, the Fluidigm logo, 48.Atlas, Dynamic Array and Juno are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2019 Fluidigm Corporation. All rights reserved. 08/2019