

# New Staining Protocols

A new protocol for staining on Helios, CyTOF 2-to-Helios Upgrade, CyTOF and CyTOF 2

## Introduction

A change to the staining protocols for routine use in all mass cytometry instruments has been developed in order to improve cell stability and staining quality. These protocols enable higher intensities and tighter CVs for signals from metal-labeled cells, improving data quality and maximizing results from precious samples. Internal testing and FAS and customer experience support improved cell stability when cells are treated with fresh fixative after antibody staining.

## Staining Protocol Updates

Updates to the staining protocols on Helios™, Helios Upgrade, CyTOF® and CyTOF 2 consist of the following:


- Addition of a fixation step using fresh fixative (1.6% FA) to sample preparation protocols prior to iridium staining in Fix/Perm, for all instruments.
- Fixative should be opened and prepared the same day it is used.

We recommend Pierce™ 16% Formaldehyde, ampule-sealed, methanol-free in 1 mL aliquots.



## Sample Fixation Prior to Iridium Staining

The addition of fresh 1.6% formaldehyde fix prior to staining with iridium staining in Fix/Perm buffer enhances the quality of your data with improved cell integrity prior to data acquisition on Helios. See the Maxpar® Cell Surface Staining Protocol with Fresh Fix (PN 400276)



**For technical support visit [techsupport.fluidigm.com](https://techsupport.fluidigm.com). | For general support visit [fluidigm.com/support](https://fluidigm.com/support).**

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [support.northamerica@fluidigm.com](mailto:support.northamerica@fluidigm.com) Latin America +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
Europe/Middle East/Africa/Russia +44 1223 859941 | [support.europe@fluidigm.com](mailto:support.europe@fluidigm.com) China (excluding Hong Kong) +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
Japan +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) All other Asian countries/India/Australia +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

**For Research Use Only. Not for use in diagnostic procedures.**

Information in this publication is subject to change without notice. **Safety data sheet information:** [fluidigm.com/sds](https://fluidigm.com/sds). **Patent and license information:** [fluidigm.com/legalnotices](https://fluidigm.com/legalnotices). **EU's WEEE directive information:** [fluidigm.com/compliance](https://fluidigm.com/compliance). Fluidigm, the Fluidigm logo, Cell-ID, EQ, Helios, and Imaging Mass Cytometry, IMC, Maxpar, and MCD are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2018 Fluidigm Corporation. All rights reserved. 04/2018