

A Method for Detecting Protein Expression in Single Cells Using the C₁TM Single-Cell Auto Prep System

Camila Egidio, Aik Ooi, Ilona Holcomb, Dave Ruff, Jing Wang, Ramesh Ramakrishnan
Fluidigm Corporation, South San Francisco, CA USA 94080

Introduction

Recent improvements in microfluidics and biochemistry have enabled single-cell molecular analysis, providing new insight into the heterogeneity of cell populations. The C₁TM Single-Cell Auto Prep System is an automated platform that streamlines the isolation and processing of 96 individual, live cells for RNA and DNA analysis. Single-cell protein profiling is a direct complement to genomic analysis as it provides additional insights into key molecular mechanisms and system biology. To enable this, we adapted a highly multiplexed protein detection method (Proseek Multiplex Oncology I^{96x96}, Olink Bioscience) based on the Proximity Extension Assay technology (PEA) for use on the C₁TM Single-Cell Auto Prep System.

Overview of the C₁TM Single-Cell Auto Prep System for Protein Detection

We have used the C₁TM Single-Cell Auto Prep System in combination with the Proximity Extension Assay technology (PEA, Figure 1A) to develop a workflow for the automated analysis of the protein expression of single cells (Figure 1B-D). The method developed is based on the use of a PEA probe panel targeting 92 different proteins and of those 66 correspond to intracellular proteins that can be detected in single cells (Figure 1E).




Figure 1

A: Schematic representation of the PEA method. Detection of the amplified reporter oligonucleotide is done by qPCR on the BioMarkTM System. Cycle threshold of the amplified reporter oligo reflects target protein abundance during the incubation step. **B:** The workflow developed for automated protein detection in single cell uses the C₁TM Single-Cell Auto Prep System that is composed of a controller instrument and integrated fluidic circuits (IFC). **C:** The C₁TM IFC architecture is shown in details, containing 96 individual capture sites and dedicated nano-chambers for downstream reactions. **D:** Representation of the system of independent chambers and valves connected to the 4.5 nL single-cell capture site in the C₁TM IFC. Each one of the 96 capture sites has its own dedicated system of chambers and valves, allowing all PEA steps to take place in a single run for 96 single cells in parallel. **C:** List of protein targets for the PEA probe panel contained in the Proseek Multiplex Oncology I^{96x96} kit used. Of the 92 protein targets, 25 (around 30%) are strictly secreted and not expected to generate signal when performing single cell analysis. **D:** Single-cell-to-result turnaround time.

Characteristic Protein Expression Signatures Identified




Figure 3

A total of 401 single cells were analyzed (represented in columns in the panel above) in eight independent C₁TM PEA experiments for each of the four human cell lines MDA-MB-231 (n=54), CRL-7163 (n=83), HL60 (n=117), and K562 (n=147) (ATCC). Protein targets are represented in horizontal lines in the panel above. Across the two experiments run for each cell line, 41, 31, 24, and 56 protein targets were detected as expressed in at least one single cell, respectively. Protein targets are considered expressed if $\Delta C_T = \text{Sample } C_T - (\text{Avg. Background } C_T - 2 \times \text{St. Dev. Background}) < -0.4$. The figure shows targets detected as expressed in a minimum of 10% of all single cells within each cell line analyzed. Of the 20 targets shown in the figure above, seven stand out as having somewhat specific expression levels in the following cell lines: Tissue Factor and IL-1ra in MDA-MB-231; Myeloperoxidase in HL60; CD69 and Cathepsin D in K562; MCP-1 and Osteoprotegerin in CRL-7163. Expression in specific cell lines and corresponding specific function was validated by literature analysis (References).

Most Targets Detected in Single Cells are Consistently Detected Across Experiments




Figure 4


The graphs above show targets detected in specific cell lines tested across two independent C₁TM PEA experiments. Since some level of variability of protein expression is expected at single-cell level, a more stringent criteria was used to select top targets expressed in the cell lines to evaluate experimental reproducibility: targets expressed in at least 10% of all single cells within at least one experiment with $\Delta C_T = \text{Sample } C_T - (\text{Avg. Background } C_T - 2 \times \text{St. Dev. Background}) < -0.4$ are shown. On average, 90% of the targets shown for each cell line were consistently expressed across the two experimental replicates at similar percentages of the cell population analyzed.

Results

Protein Detection in Single Cells is Consistent Across Chip and Plate Experiments

Figure 5

Results from PEA on plate-sorted cells were compared to results obtained from two independent C₁TM PEA experiments on single HL60 cells. In general, results obtained from plate PEA on sorted cells confirmed results obtained by C₁TM PEA, with the exception of Tissue Factor. However, plate PEA signal for this specific target does not increase as expected when 10 and 50 cells are tested, suggesting that the high background signal of plate PEA could be affecting expression level results for this method.



Flow Cytometry and Immunofluorescence Results are Consistent with C₁TM PEA Results




Figure 6

A: C₁TM PEA results for two specific targets were validated on HL60 and K562 cells using orthogonal methods: EpCAM (low and high expression, respectively) and EMMPRIN (high expression in both cell types) antibodies conjugated with fluorescent dyes were used to evaluate expression levels of populations of cells with flow cytometry (Flow) and for on-chip immunofluorescence (IF) on single cells prior to C₁TM PEA. Flow and IF results were highly concordant with PEA results and some of the expression rate differences observed can be explained by different antibodies used across the methods and different population of cells tested (flow vs. PEA and IF). **B:** The diagram shows a heat map of the protein expression results for C₁TM PEA and IF for EpCAM (red indicates high expression). As expected, K562 cells have high EpCAM expression confirmed by PEA and IF and HL60 cells have high MPO expression levels confirmed by PEA. Two cells out of 38 analyzed with IF and PEA had results different than expected, presenting both EpCAM expression (IF and PEA) and MPO (PEA). For one of those cells it has been confirmed that two instead of one cell had been captured in the C₁TM IFC chamber (panel C).

Conclusion

- We have developed methodology for automated protein detection from single cells on the C₁TM Single-Cell Auto Prep System, with the ability to simultaneously process up to 96 single cells.
- The method is sensitive enough to detect expression levels from single cells and is a promising technique to use in combination with DNA and RNA profiling from single cells for further system biology studies. It is also consistent with other studies that target gene expression (References).
- The PEA probe panel from the Proseek Multiplex Oncology I^{96x96} kit, which targets 92 potential cancer-related targets, has been successfully used in profiling single cells derived from both cancer and normal tissue, grouping 98% of all cells analyzed (n=401).

References

- Fang et al. "The role of a new CD44st in increasing the invasion capability of the human breast cancer cell line MCF-7". BMC Cancer 11:290 (2011).
- Van Lint et al. "Chemokine and cytokine processing by matrix metalloproteinases and its effect on leukocyte migration and inflammation". J Leuk Bio 82:6 (2007); 1375-1381.
- Yao et al. "Important Functional Roles of Basigin in Thymocyte Development and T cell Activation". Int J Biol Sci 10(1):43-52 (2014).
- O'Donovan et al. "Caspase 3 in Breast Cancer". Clin Cancer Res 9:738 (2003).
- Doerfler et al. "Caspase Enzyme Activity Is Not Essential for Apoptosis During Thymocyte Development". J Immunol 164:8 (2000) 4071-4079.
- Munz et al. "The carcinoma-associated antigen EpCAM upregulates c-myc and induces cell proliferation". Oncogene 23 (34): 5748-58 (2004).
- Versteeg et al. "Tissue Factor and Cancer Metastasis: The Role of Intracellular and Extracellular Signaling Pathways". Mol Med 10(1-6): 6-11 (2004 Jan-Jun).
- Murao et al. "Myeloperoxidase: a myeloid cell nuclear antigen with DNA-binding properties". PNAS 85(4): 1232-1236 (1988).
- Hantschel et al. "The chemokine interleukin-8 and the surface activation protein CD69 are markers for Bcr-Abl activity in chronic myeloid leukemia". Mol Oncol 2(3):272-81 (2008 Oct).
- Lkholder et al. "Cathepsin D released by lactating rat mammary epithelial cells is involved in prolactin cleavage under physiological conditions". J Cell Science 117 (Pt 21): 5155-5164.
- Burn et al. "Monocyte chemoattractant protein-1 gene is expressed in activated neutrophils and retinoic acid-induced human myeloid cell lines". Blood 84:8 2776-2783 (1994).
- Fisher et al. "Osteoprotegerin over-expression by breast cancer cells enhances orthotopic and osseous tumor growth, and contrasts with that delivered therapeutically". Cancer Research 66, 3620-3628.