

Callisto Adherent Cell Culture Reagent Kit, 5 IFCs

Kit Part Number: 100-8682

Intended Use

The purpose of this kit is to provide the necessary materials to support the Callisto™ system.

Contents

PN	Reagent	Quantity	Volume/Unit (mL)	Uses/Unit	Storage
100-8505	Callisto Blocking Reagent	5 tubes	1	5 IFCs	4 °C
100-8506	Cell Rinsing Reagent	5 tubes	1.6	5 IFCs	
100-8507	Callisto Suspension Reagent	5 tubes	0.5	5 IFCs	
100-8508	Callisto Pre-Harvest Reagent	5 tubes	1	5 IFCs	
100-8504	Callisto Control Line Fluid	1 bottle	10	25 uses	
100-6250	Callisto Actuation Fluid	2 tubes	1.5	5 IFCs	

Handling

Before using reagents, refer to the Callisto System Getting Started Guide (PN 100-7599).

Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

For technical support visit fluidigm.com/support

North America +1 650 266 6100 | Toll-free: +1 866 358 4354 in the US | techsupport@fluidigm.com
Europe +33 1 60 92 42 40 | techsupporteurope@fluidigm.com
China (excluding Hong Kong) +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com
All other Asian countries +1 650 266 6100 | techsupportasia@fluidigm.com
Central and South America +1 650 266 6100 | techsupportlatam@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Safety data sheet information:** fluidigm.com/sds. **Patent and license information:** fluidigm.com/legalnotices. **EU's WEEE directive information:** fluidigm.com/compliance. Fluidigm, the Fluidigm logo, and Callisto are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2015 Fluidigm Corporation. All rights reserved. 07/2015