

# Advanta NGS Library Prep Reagent Bundle— LP 48.48, 10 IFCs








Part Number: 101-7666

## Purpose


The Advanta™ NGS Library Prep Reagent Bundle—LP 48.48, 10 IFCs provides the necessary reagents to support next-generation sequencing (NGS) library preparation on the LP 48.48 IFC. The kit contains sufficient reagents for 10 integrated fluidic circuits (IFCs).

## Contents

### 101-7665, Advanta NGS Library Prep Reagent Kit—LP 48.48, 10 IFCs

Part Number	Product Name	Cap Color	Unit Quantity	Volume/Unit	Storage
101-7634	TSP Sample Loading Reagent v2	 Purple	1 tube	200 µL	
101-0408	TSP Adapter Mix	 Light blue	1 tube	120 µL	
100-5941	PCR Water	 Natural	2 tubes	1.8 mL	
100-9167	DNA Dilution Reagent	 Green	1 bottle	3.0 mL	–20 °C
101-0409	TSP Assay Loading Reagent	 Light tan	2 tubes	75 µL	
101-3050	4X TSP Master Mix	 Yellow	1 tube	1.05 mL	
101-0995	TSP DNA Polymerase	 Orange	4 tubes	50 µL	

### 101-0741, Targeted DNA Seq Harvest Reagent—5 Bottles

Part Number	Product Name	Cap Color	Unit Quantity	Volume/Unit	Storage
101-0743	TSP Harvest Reagent	 Clear	5 bottles	10 mL	–20 °C

## Handling

Before using reagents, refer to the Advanta NGS Library Preparation with Juno™ Protocol (PN 101-7878) or the Advanta NGS Library Preparation with Access Array™ Protocol (PN 101-7885).

## Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to [fluidigm.com/sds](http://fluidigm.com/sds) and search for the SDS using either the product name or the part number.

**For technical support visit [fluidigm.com/support](http://fluidigm.com/support).**

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com) Latin America +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
 Europe/Middle East/Africa/Russia +44 1223 859941 | [techsupporteurope@fluidigm.com](mailto:techsupporteurope@fluidigm.com) China (excluding Hong Kong) +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
 Japan +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) All other Asian countries/India/Australia +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

**For Research Use Only. Not for use in diagnostic procedures.**

Information in this publication is subject to change without notice. **Patent and license information:** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). **Trademarks:** Fluidigm, the Fluidigm logo, Access Array, Advanta and Juno are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries.

© 2018 Fluidigm Corporation. All rights reserved. 05/2018