

Air Filter, Large

Catalog#: 101430
Package Size: 1 unit
Instrument: CyTOF

Technical Information

Description:

The Large Air filter is necessary for preventing solid particles such as dust from the air from entering the CyTOF® Mass Cytometer.

The Large Air Filter is installed on the left rear side of the instrument.



Contact Information:

Sales: sales@DVSsciences.com | Support: Support.NorthAmerica@DVSsciences.com or Support.Europe@DVSsciences.com
www.DVSsciences.com | For assistance by phone: +1-855-DVS-CYTO

This product is intended for Research Use Only. This product is not intended for therapeutic or diagnostic purposes in humans or animals.

Vendor info & description quote as needed.

List supplier of raw material antibody here only if targets/clones are unique, licensed, and/or indicated as required in our agreement.

© 2011 DVS Sciences Inc. The trademarks mentioned herein are the property of DVS Sciences Inc. or their respective owners.

DVS Sciences Inc. products are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to 6 months from date of sales for products used, handled and stored according to DVS Sciences, Inc. instructions. DVS Sciences, Inc. sole liability for the product is limited to replacement of the product or refund of the purchase price. Product may not be resold, modified for resale or used to manufacture commercial products without prior written approval from DVS Sciences Inc.

