

## Anti-Phospho-Rb[pS807/pS811]-166Er

**Catalog #:** 3166011A

**Package Size:** 50 tests

**Storage:** Store product at 4°C. Do not freeze.

**Cross Reactivity:** Rat, Mouse, Human

**Clone:** J112-906


**Isotype:** Mouse IgG2a

**Formulation:** Antibody stabilizer with 0.05% Sodium Azide

### Technical Information

**Validation:** Each lot of conjugated antibody is quality control tested by CyTOF® analysis of stained cells using the appropriate positive and negative cell staining and/or activation controls.

**Recommended Usage:** The suggested use is 1 µl for up to 3 X 10<sup>6</sup> live cells in 100 µl. It is recommended that the antibody be titrated for optimal performance for each of the desired applications.



Human T-lymphoblasts MOLT-4 cells were fixed, permeabilized, and stained with 162Dy-anti-Ki-67 (B56) and 166Er-anti-pRb [Ser807/811](J112-906).

### Description

The Retinoblastoma (Rb) gene encodes a phosphoprotein that negatively regulates the cell cycle, and this activity is critical for the classic RB-mediated tumor suppression function. Mutations in the retinoblastoma (RB) susceptibility gene lead to childhood retinal cancer, and the RB pathway is altered in 70% of human cancer types. Upon mitogenic signaling, RB protein is inhibited by hyperphosphorylation, which disrupts RB transcriptional repression complexes to allow G1 to S phase transition. RB also induces S-phase arrest, triggered by the intra S-phase checkpoint in the presence of DNA-damaging agents, such as cisplatin and VP-16. RB also participates in other cellular processes, such as terminal cell differentiation and maintenance of genetic stability. It has also been suggested that RB can exert an anti-apoptotic function, and that RB is cleaved by caspases during apoptosis induced by VP-16 and Tumor Necrosis Factor- $\alpha$  (TNF- $\alpha$ ).

### References

Bandura, D. R., et al. Mass Cytometry: Technique for Real Time Single Cell Multitarget Immunoassay Based on Inductively Coupled Plasma Time-of-Flight Mass Spectrometry. *Analytical Chemistry* 81:6813-6822, 2009.

Behbehani, G.K., et al. Single-cell mass cytometry adapted to measurements of the cell cycle. *Cytometry A* 81 (7): 552-566, 2012.

Ornatsky, O. I., et al. Highly multiparametric analysis by mass cytometry. *J Immunol Methods* 361 (1-2):1-20, 2010.

### Contact Information:

Sales: sales@DVSsciences.com | Support: support@DVSsciences.com  
www.DVSsciences.com | For assistance by phone: 855-DVS-CYTO

This product is intended for Research Use Only. This product is not intended for therapeutic or diagnostic purposes in humans or animals.

© 2011 DVS Sciences Inc. The trademarks mentioned herein are the property of DVS Sciences Inc. or their respective owners. DVS Sciences Inc. products are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to 6 months from date of sales for products used, handled and stored according to DVS Sciences, Inc. instructions. DVS Sciences, Inc. sole liability for the product is limited to replacement of the product or refund of the purchase price. Product may not be resold, modified for resale or used to manufacture commercial products without prior written approval from DVS Sciences Inc.

