

# Flex Six SNP Type Genotyping Reagent Kit—5 IFCs

**Kit Part Number: 100-7355**

## Purpose

The Flex Six™ SNP Type Genotyping Reagent Kit provides necessary reagents to support genotyping using Flex Six Genotyping IFCs with SNP Type™ assays. The kit contains sufficient reagents for 5 integrated fluidic circuits (IFCs).

## Contents

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
100-5359	2X Assay Loading Reagent	● Yellow	1.0 mL/tube	2 tubes	
100-7445	20X SNP Type Sample Loading Reagent	○ White	68 µL/tube	2 tubes	-20 °C
100-3402	60X SNP Type Reagent	● Amber	70 µL/tube	1 tube	
101-5941	PCR Water	○ Natural	1.8 mL/tube	1 tube	

## Handling

Before using this kit, see the SNP Genotyping Analysis User Guide (PN 68000098).

## Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to [fluidigm.com/sds](http://fluidigm.com/sds) and search for the SDS using either the product name or the part number.

### For technical support visit [techsupport.fluidigm.com](http://techsupport.fluidigm.com).

**North America** +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com) **Latin America** +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
**Europe/Middle East/Africa/Russia** +44 1223 598100 | [techsupporteurope@fluidigm.com](mailto:techsupporteurope@fluidigm.com) **China (excluding Hong Kong)** +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
**Japan** +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) **All other Asian countries/India/Australia** +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

### For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). Fluidigm, the Fluidigm logo, Flex Six, and SNP Type are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2018 Fluidigm Corporation. All rights reserved. 12/2018