

4X Assay Loading Reagent—5 Tubes

Kit Part Number: 102-0135

Purpose

The 4X Assay Loading Reagent kit provides a reagent to support real-time PCR experiments using Dynamic Array™ integrated fluidic circuits (IFCs).

Contents

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
102-0114	4X Assay Loading Reagent	● Yellow	120 µL/tube	5 tubes	-20 °C

Handling

Before using the reagent, see the Real-Time PCR for Viral RNA Detection Quick Reference (FLDM-00102).

Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

For technical support visit techsupport.fluidigm.com.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com **Latin America** +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 598100 | techsupporteurope@fluidigm.com **China** (excluding Hong Kong) +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com **All other Asian countries/India/Australia** +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legal/notices. **Trademarks:** Fluidigm, the Fluidigm logo, and Dynamic Array are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2020 Fluidigm Corporation. All rights reserved. 04/2020