

# Advanta PCR MM—4 Tubes

**Kit Part Number: 102-0876**

## Purpose

The Advanta™ PCR MM kit provides a reagent to support the Advanta RT-PCR Kit.

## Contents

Part Number	Component	Cap Color	Volume/Unit	Unit Quantity	Storage
102-0420	Advanta PCR MM	● Red	1.2 mL/tube	4 tubes	-20 °C

## Handling

Before using the assays, see the Pathogen Detection Using the Advanta RT-PCR Kit Protocol (FLDM-00191).

## Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to [fluidigm.com/sds](https://fluidigm.com/sds) and search for the SDS using either the product name or the part number.

### For technical support visit [techsupport.fluidigm.com](https://techsupport.fluidigm.com).

**North America** +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com) **Latin America** +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
**Europe/Middle East/Africa/Russia** +33 1 60 92 42 40 | [eu.support@fluidigm.com](mailto:eu.support@fluidigm.com) **China** (excluding Hong Kong) +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
**Japan** +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) **All other Asian countries/India/Australia** +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

### For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** [fluidigm.com/legal/notices](https://fluidigm.com/legal/notices). **Trademarks:** Fluidigm, the Fluidigm logo, and Advanta are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2020 Fluidigm Corporation. All rights reserved. 11/2020